



YOUR
DISPOSABLE
EXCELLENCE

30
YEARS

Protect, Perform, Personalize!

Dispo Roll: disposable medical sheets designed for you



Dispo Roll is
versatile and
customizable
in every detail:



COLOR




DIMENSION



GRAMMAGE

Dispo Roll finds
extensive use in
these and many
other areas





YOUR
DISPOSABLE
EXCELLENCE

30
YEARS

Features of Dispo Roll Disposable Medical Sheets

- Dispo Roll disposable medical sheets are composed of an absorbent layer of cellulose wadding, a natural and soft material, and a layer of waterproof polyethylene.
- They have excellent breathability, which helps keep the patient's skin dry and fresh.




Technical information for Dispo Roll available in the catalog:

Color: blue or white, Height: 50, 60 and other on request

Length: from 50 to 100 mt, Perforation: every 38 cm

Grammage of cellulose: 23 or 26 g/m²

PE film: available in HDPE 7 µm or in LDPE 12 µm



For more
Information



DISPOTECH s.r.l.

Via Al Piano, 29 – 23020 Gordona (So)

+39 0343 36711 – info@dispotech.it